



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 02/05/MO	Status 02/05/MO
Substance code:		Page : 1 (3)
1. Identification of the substance / preparation and company		
Product details		
Trade Name: C&B Temp Liquid		
Supplier details :		
PSP Dental Company Limited 3-5 Dylan Road Belvedere Kent DA17 5QS		Information provided by: PSP Dental Company Limited
		Emergency telephone number: 020 8311 7337 Emergency fax number: 020 8310 0920
2. Composition/information on ingredients		
Chemical characterization		
Hazardous Components		
Ingredient	%	CAS No. S Phrase R Phrase
Iso Butyl Methacrylate	48.55	97-86-9 S 29 R 38
Ethyl Methacrylate	43.70	97-63-2 S 25 R 37
Others	7.79	- -
3. Hazards identification		
Unlikely to cause harmful effects under normal conditions of handling and use. Practice good hygiene. Wear gloves, overalls.		
4. First aid measures		
After contact with skin: Mild irritant, prolonged contact may cause defatting and dermatitis. Wash skin first with water and then, soap and water. If in doubt seek medical help.		
After contact with eyes: Mild irritant, wash open eye with of clean water for 10 minutes. Seek medical help.		
After ingestion: Seek medical advice.		
After inhalation: Remove to fresh air. If breathing is irregular or stopped, give artificial respiration. Give nothing by mouth. Seek medical help.		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 2 (3)
5. Fire-fighting measures: Water spray, foam, dry powder or carbon dioxide. Self contained breathing apparatus.		
6. Accidental release measures: Environmental precautions Prevent entry into drains and water courses. Avoid creating dusts. Methods for cleaning up / taking up: Wear approved respirator suitable for chemical particle concentrations encountered and impervious gloves (for bulk spillages). Sweep spill into containers for disposal. Wash contaminated area, wash affected skin areas with soap and water.		
7. Handling and storage Handling: Good General Ventilation Follow instructions enclosed in package and use only for purposes directed. Avoid contact with eyes. Avoid creating dusts. Storage: The container should be stored in normal dry place with lid firmly closed, since material is moisture sensitive.		
8. Exposure controls / personal protection Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work. Respiratory protection: Type approved RPE for organic vapour if required. Hand protection: rubber gloves. Eye protection: Goggles. Body protection: Overalls and boots		
9. Physical and chemical properties Appearance Form : Liquid Colour : Clear and colourless liquid Odour : Ester like Data relevant to safety Flash point : 43EC Ignition temperature : N/A Explosive Limits : 1.8 vol % up to saturation Boiling Point / Range : 1013 @ + 155EC Melting Point Range: Density : 0.888 (20EC) Solubility in water: Negligible Vapour Pressure: 4 @ 20EC		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 3 (3)
10. Stability and reactivity Thermal decomposition / Conditions to avoid Stable in normal conditions Hazardous reactions / Substances to avoid Polymerization may be initiated by contamination with peroxides, azo compounds, heavy metal ions, test.amines. Polymerization is also induced by light.		
11. Toxicological information There is no Toxicological data available for this product. There is no known health risk if used properly. Prolonged skin and eye contact should be avoided. Remarks		
12. Ecological information Further ecological information Remarks		
13. Disposal considerations Product Substance via authorized waste disposal contractor to an approved waste disposal site, observing all local and national regulations.		
14. Transport information UN No: 2277 IMO: N/A IATA: 2277 IMDG Class: Cl 3.2 Packaging Group: N/A Packaging Group: N/A Correct Technical Name: ADR/RID: Cl. 3.36 Not regulated for transport		
15. Regulatory information Hazard symbols: R phrases: R 37 S phrases: S25: Avoid contact with eyes Other Regulations:		
16. Other information This Information is given in good faith, being based on the latest knowledge available to P.S.P.Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P.Dental Co. Ltd prior to use for further advice.		