



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 05/05/MO	Status 05/05/MO		
Substance code:		Page : 1 (3)		
1. Identification of the substance / preparation and company				
Product details				
Trade Name: C & B Temp Powder				
Supplier details :				
PSP Dental Company Limited		Information provided by: PSP Dental Company Limited		
3-5 Dylan Road				
Belvedere				
Kent DA17 5QS				
Emergency telephone number: 020 8311 7337				
Emergency fax number: 020 8310 0920				
2. Composition/information on ingredients				
Chemical characterization				
Hazardous Components				
Ingredient	%	CAS No.	S Phrase	R Phrase
Ethyl Methacrylate	98.5	000094360	S24	R 43
Others	1.5	-	-	-
3. Hazards identification				
Unlikely to cause harmful effects under normal conditions of handling and use.				
Practice good hygiene. Wear gloves, overalls.				
4. First aid measures				
After contact with skin:				
Mild irritant, prolonged contact may cause defatting and dermatitis. Wash skin first with water and then, soap and water. If in doubt seek medical help.				
After contact with eyes:				
Mild irritant, wash open eye with of clean water for 10 minutes. Seek medical help.				
After ingestion:				
Seek medical advice.				
After inhalation:				
Remove to fresh air. If breathing is irregular or stopped, give artificial respiration. Give nothing by mouth. Seek medical help.				



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 05/05/MO	Status	05/05/MO
Substance code:		Page :	2 (3)
5. Fire-fighting measures: Normal fire fighting procedures may be used			
6. Accidental release measures: Environmental precautions Prevent entry into drains and water courses. Avoid creating dusts. Methods for cleaning up / taking up: Wear approved respirator suitable for chemical dust concentrations encountered and impervious gloves (for bulk spillages). Sweep spill into containers for disposal. Wash contaminated area, wash affected skin areas with soap and water.			
7. Handling and storage Handling: Good General Ventilation Follow instructions enclosed in package and use only for purposes directed. Avoid contact with eyes. Avoid creating dusts. Storage: The container should be stored in normal dry place with lid firmly closed, since material is moisture sensitive.			
8. Exposure controls / personal protection Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work. Respiratory protection: Type approved RPE for dusts if required. Hand protection: rubber gloves. Eye protection: Goggles. Body protection: Overalls and boots			
9. Physical and chemical properties Appearance Form : Powder Colour : Light Yellow Odour : Distinct Data relevant to safety Flash point : Melting Point Range: 150 - 230 EC Ignition temperature : N/A Density : 550 - 750 Kg / m3 Explosive Limits : N/A Solubility in water: Forms paste with water. Boiling Point / Range : N/A Vapour Pressure: N/A			



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 3 (3)
10. Stability and reactivity	Thermal decomposition / Conditions to avoid Stable in normal conditions	
	Hazardous reactions / Substances to avoid Unstable above 98EC gives Acrylate Vapour	
11. Toxicological information	There is no Toxicological data available for this product. There is no known health risk if used properly. Prolonged skin and eye contact should be avoided. Remarks	
12. Ecological information	Further ecological information Remarks	
13. Disposal considerations	Product Substance via authorized waste disposal contractor to an approved waste disposal site, observing all local and national regulations.	
14. Transport information	Not classified as dangerous for transport	
	UN No:	IMDG Class:
	IMO:	Packaging Group:
	IATA:	Packaging Group:
	Correct Technical Name:	
	ADR/RID:	
	Not regulated for transport	
15. Regulatory information	Hazard symbols:	
	R phrases: R 43: May cause sensitisation by skin contact	
	S phrases: S25: Avoid contact with eyes	
	Other Regulations: This data was prepared in accordance with Directive 2001/58/EC	
16. Other information	This Information is given in good faith, being based on the latest knowledge available to P.S.P.Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P.Dental Co. Ltd prior to use for further advice.	