



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/05/YB	Page :	1 (3)
Substance code:			
1. Identification of the substance / preparation and company			
Product details			
Trade Name: Ceram II VLC Liquid			
Supplier details :			
PSP Dental Company Limited		Information provided by:	
3-5 Dylan Road		PSP Dental Company Limited	
Belvedere			
Kent DA17 5QS			
		Emergency telephone number: 020-8311-7337	
		Emergency fax number: 020 8310 0920	
2. Composition/information on ingredients			
Chemical characterization			
Hazardous Components			
Ingredient	%	CAS No.	S Phrase R Phrase
HEMA		10-30 868-77-9	
3. Hazards identification			
Unlikely to cause harmful effects under normal conditions of handling and use.			
Practice good hygiene. Wear gloves, overalls.			
4. First aid measures			
After contact with skin:			
Slight irritation (redness) upon repeated or prolonged contact with skin.			
After contact with eyes:			
Irritant, wash open eye with clean water for 10 minutes. Seek medical help.			
After ingestion:			
Seek medical advice immediately, do not make the casualty vomit. If conscious, rinse mouth with water and give water to drink. Repeat every 10 minutes.			
After inhalation:			
N/A			



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/05/YB	Page :	2 (3)
Substance code:			
5. Fire-fighting measures: Foam, Dry Powder, waste spray/fog. As with all fires involving chemicals, chemical-resistant protective suit, gloves, boots and breathing apparatus should be worn.			
6. Accidental release measures: Environmental precautions Prevent entry into drains and water courses. Methods for cleaning up / taking up: Dry mop, then clean with water. Use protective clothing and gloves.			
7. Handling and storage Handling: Good General Ventilation Follow instructions enclosed in package and use only for purposes directed. Avoid contact with eyes. Storage: The container should be stored in normal dry place with lid firmly closed, since material is light sensitive			
8. Exposure controls / personal protection Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work. Respiratory protection: N/A Hand protection: rubber gloves. Eye protection: Goggles. Body protection: Overalls and boots			
9. Physical and chemical properties Appearance Form : Liquid Colour : Light Yellow Odour : Acrylic Data relevant to safety Flash point : > 100EC Melting Point Range: N/A Ignition temperature : N/A Density : N/A Explosive Limits : N/A Solubility in water: Miscible Boiling Point / Range : > 200EC Vapour Pressure: N/A			



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/05/YB	Page :	3 (3)
Substance code:			
10. Stability and reactivity	Thermal decomposition / Conditions to avoid Stable in normal conditions Incompatibility: None Hazardous reactions / Substances to avoid Reactive materials, e.g: metals, peroxides etc. Temperature greater than 40EC and light		
11. Toxicological information	LD ₅₀ (Rat) > 200mg/kg There is no known health risk if used properly. Prolonged skin and eye contact should be avoided.		
12. Ecological information	Further ecological information Remarks: N/A		
13. Disposal considerations	Product Substance via authorized waste disposal contractor to an approved waste disposal site, observing all local and national regulations.		
14. Transport information	UN No: IMDG Class: IMO: Packaging Group: IATA: Packaging Group: Correct Technical Name: ADR/RID: Not regulated for transport		
15. Regulatory information	Hazard symbols: R phrases: S phrases: S25: Avoid contact with eyes Other Regulations:		
16. Other information	This Information is given in good faith, being based on the latest knowledge available to P.S.P.Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P.Dental Co. Ltd prior to use for further advice.		