



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 09/98/MW	Status 09/98/MW
Substance code:		Page : 1 (3)
1. Identification of the substance / preparation and company		
Product details : Glass Ionomer Core Restorative Silver Powder (Aqua Set)		
Trade Name: Ceramcore B Silver		
Supplier details :		
PSP Dental Company Limited 3-5 Dylan Road Belvedere Kent DA17 5QS		Information provided by: PSP Dental Company Limited
		Emergency telephone number: 020-8311-7337 Emergency fax number: 020 8310 0920
2. Composition/information on ingredients		
Chemical characterization		
Ingredient		
%	CAS No.	S Phrase R Phrase
	Verified mixture of Calcium/Soduin, Fluoride	71.0
	Aluminium Oxides/Phosphates and silicates	18.0
	Polyacrylic Acid	8.0
	Silver Alloy	2.0
	Retarder/Pigment	Trace
3. Hazards identification		
Unlikely to cause harmful effects under normal conditions of handling and use. Practice good hygiene. Wear gloves, overalls.		
4. First aid measures		
After contact with skin: Wash skin thoroughly with soap and water, remove and wash clothing before use.		
After contact with eyes: Mild irritant, wash open eye with plenty of clean water for 15 minutes. Seek medical help if irritation persists.		
After ingestion: Accidental ingestion of hardened Ceramcore Silver filling is not known to cause any health problems. Ingestion of a large quantity of Ceramcore Silver powder will cause gastro irritation and immediate medical care should be sought.		
After inhalation: Remove to fresh air. If breathing is irregular or stopped, give artificial respiration. Give nothing by mouth. Seek medical help if discomfort persists.		

