



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 02/06/YB	Status
Substance code:		Page : 1 (3)
1. Identification of the substance / preparation and company		
Product details		
Trade Name: Ceramlin B (Aqua Set) Powder		
Supplier details :		
PSP Dental Company Limited		Information provided by: PSP Dental Company Limited
3-5 Dylan Road		
Belvedere		
Kent DA17 5QS		
Emergency telephone number: 020 8311 7337		
Emergency fax number: 020 8310 0920		
2. Composition/information on ingredients		
Chemical characterization		
Hazardous Components		
Ingredient	%	CAS No. S Phrase R Phrase
Polyacrylic Acid	21.0	
Alumino Silicate Fluoride Barium Glass	78.95	
Iron Oxide Pigment (approved)	0.035	
3. Hazards identification		
Unlikely to cause harmful effects under normal conditions of handling and use.		
Practice good hygiene. Wear gloves, overalls.		
4. First aid measures		
After contact with skin:		
Mild irritant. Wash skin first with water and then, soap and water. If in doubt seek medical help.		
After contact with eyes:		
Mild irritant, wash open eye with clean water for 15 minutes. Seek medical help.		
After ingestion: Seek medical advice.		
After inhalation:		
Remove to fresh air. If breathing is irregular or stopped, give artificial respiration. Give nothing by mouth. Seek medical help.		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 02/06/YB	Status
Substance code:		Page : 2 (3)
5. Fire-fighting measures: Suitable extinguishing media. Material is non combustible. Use extinguishing method and media suitable to surrounding fire.		
6. Accidental release measures: Environmental precautions Prevent entry into drains and water courses. Avoid creating dusts. Methods for cleaning up / taking up: Wear approved respirator suitable for chemical dust concentrations encountered and impervious gloves (for bulk spillages). Sweep spill into containers for disposal. Wash contaminated area, wash affected skin areas with soap and water.		
7. Handling and storage Handling: Good General Ventilation Follow instructions enclosed in package and use only for purposes directed. Avoid contact with eyes. Avoid creating dusts. Storage: The container should be stored in normal dry place with lid firmly closed, since material is moisture sensitive.		
8. Exposure controls / personal protection Hygiene measures Remove soiled clothing immediately. Wash hands before breaks and after work. Respiratory protection: Type approved RPE for dusts if required. Hand protection: rubber gloves. Eye protection: Goggles. Body protection: Overalls and boots		
9. Physical and chemical properties Appearance Form : Powder (Fine) Colour b : Light Yellow Odour : Typical Data relevant to safety Flash point : N/A Melting Point Range: N/A Ignition temperature N/A Density : N/A Explosive Limits : N/A Solubility in water: Forms paste with water. Boiling Point / Range : N/A Vapour Pressure:		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 02/06/YB	Status
Substance code:		Page : 3 (3)
10. Stability and reactivity	Thermal decomposition / Conditions to avoid Stable in normal conditions Hazardous reactions / Substances to avoid Unstable above 98EC gives Acrylate Vapour	
11. Toxicological information	There is no Toxicological data available for this product. There is no known health risk if used properly. Prolonged skin and eye contact should be avoided. Remarks	
12. Ecological information	Further ecological information Remarks	
13. Disposal considerations	Product Substance via authorized waste disposal contractor to an approved waste disposal site, observing all local and national regulations.	
14. Transport information	UN No: N/A IMO: N/A IATA: N/A IMDG Class: N/A Packaging Group: N/A Packaging Group: N/A Correct Technical Name: ADR/RID: Not regulated for transport	
15. Regulatory information	Hazard symbols: R phrases: S phrases: S25: Avoid contact with eyes Other Regulations:	
16. Other information	This Information is given in good faith, being based on the latest knowledge available to P.S.P.Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P.Dental Co. Ltd prior to use for further advice.	