



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 09/98/MW	Status 09/98/MW
Substance code:		Page : 1 (3)
1. Identification of the substance / preparation and company		
Product details		
Trade Name: Comf Dent Liquid		
Supplier details :		
PSP Dental Company Limited 3-5 Dylan Road Belvedere Kent DA17 5QS		Information provided by: PSP Dental Company Limited
		Emergency telephone number: 020 8311 7337 Emergency fax number: 020 8310 0920
2. Composition/information on ingredients		
Chemical characterization		
Hazardous Components		
Ingredient	%	CAS No. S Phrase R Phrase
Butyl Phthalyl Butyl (Glycolate)	93.0	85-70-1
Minor constituents	6.5	
3. Hazards identification		
Unlikely to cause harmful effects under normal conditions of handling and use. Practice good hygiene. Wear gloves, overalls.		
4. First aid measures		
After contact with skin:		
Mild irritant, prolonged contact may cause defatting and dermatitis. Wash skin first with water and then, soap and water. If in doubt seek medical help.		
After contact with eyes:		
Mild irritant, wash open eye with of clean water for 15 minutes. Seek medical help.		
After ingestion:		
Seek medical advice.		
After inhalation:		
Remove to fresh air. If breathing is irregular or stopped, give artificial respiration. Give nothing by mouth. Seek medical help.		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 2 (3)
5. Fire-fighting measures: Klear self-contained breathing apparatus and protective clothing. Normal fire fighting procedures may be used.		
6. Accidental release measures: Environmental precautions Prevent entry into drains and water courses. Avoid creating dusts. Methods for cleaning up / taking up: Contain spilled liquid with sand or vermiculite and place in chemical waste container		
7. Handling and storage Handling: Good General Ventilation Follow instructions enclosed in package and use only for purposes directed. Avoid contact with eyes, skin and clothing . Use with adequate ventilation. Wash thoroughly after handling. Storage: The container should be stored in normal dry place with lid firmly closed, since material is moisture sensitive.		
8. Exposure controls / personal protection Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work. Respiratory protection: The use of a NIOSH-approved organic vapour respirator is recommended Hand protection: rubber gloves. Eye protection: Goggles. Body protection: Overalls and boots		
9. Physical and chemical properties Appearance Form : Clear, oily liquid Colour : Clear Odour : Slight characteristic Data relevant to safety Flash point : Melting Point Range: less than 35EC Ignition temperature : N/A Density : N/A Explosive Limits : N/A Solubility in water: <input type="checkbox"/> 0.1 % Boiling Point / Range : N/A Vapour Pressure: 3.9 mm Hg @ 200EC		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 3 (3)
10. Stability and reactivity	Thermal decomposition / Conditions to avoid Stable in normal conditions	
	Hazardous reactions / Substances to avoid As with other organic materials, combustion will produce carbon dioxide (CO ₂) and Carbon monoxide (CO)	
11. Toxicological information	There is no Toxicological data available for this product. There is no known health risk if used properly. Prolonged skin and eye contact should be avoided. Remarks	
12. Ecological information	Further ecological information Remarks : Is expected to rapidly hydrolyze in environment media	
13. Disposal considerations	Product Substance via authorized waste disposal contractor to an approved waste disposal site, observing all local and national regulations.	
14. Transport information	UN No: IMO: N/A IATA: Correct Technical Name: ADR/RID: Not regulated for transport	
	IMDG Class:	Packaging Group: Packaging Group:
15. Regulatory information	Hazard symbols: R phrases: S phrases: S25: Avoid contact with eyes Other Regulations:	
16. Other information	This Information is given in good faith, being based on the latest knowledge available to P.S.P.Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P.Dental Co. Ltd prior to use for further advice.	