



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 05/05/MO	Status 05/05/MO		
Substance code:		Page : 1 (3)		
1. Identification of the substance / preparation and company				
Product details				
Trade Name: Comf- Dent Powder				
Supplier details :				
PSP Dental Company Limited 3-5 Dylan Road Belvedere Kent DA17 5QS		Information provided by: PSP Dental Company Limited		
		Emergency telephone number: 020 8311 7337 Emergency fax number: 020 8310 0920		
2. Composition/information on ingredients				
Chemical characterization				
Hazardous Components				
Ingredient	%	CAS No.	S Phrase	R Phrase
Ethyl Methacrylate	99.5	0000 94 360	S 24	R 43
Others	0.47	-	-	-
3. Hazards identification				
Unlikely to cause harmful effects under normal conditions of handling and use. Practice good hygiene. Wear gloves, overalls.				
4. First aid measures				
After contact with skin:				
Mild irritant, prolonged contact may cause defatting and dermatitis. Wash skin first with water and then, soap and water. If in doubt seek medical help.				
After contact with eyes:				
Mild irritant, wash open eye with of clean water for 10 minutes. Seek medical help.				
After ingestion:				
Seek medical advice.				
After inhalation:				
Remove to fresh air. If breathing is irregular or stopped, give artificial respiration. Give nothing by mouth. Seek medical help.				



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 05/05/MO	Status 05/05/MO
Substance code:		Page : 2 (3)
5. Fire-fighting measures: Normal fire fighting procedures may be used		
6. Accidental release measures: Environmental precautions Prevent entry into drains and water courses. Avoid creating dusts. Methods for cleaning up / taking up: Wear approved respirator suitable for chemical dust concentrations encountered and impervious gloves (for bulk spillages). Sweep spill into containers for disposal. Wash contaminated area, wash affected skin areas with soap and water.		
7. Handling and storage Handling: Good General Ventilation Follow instructions enclosed in package and use only for purposes directed. Avoid contact with eyes. Avoid creating dusts. Storage: The container should be stored in normal dry place with lid firmly closed, since material is moisture sensitive.		
8. Exposure controls / personal protection Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work. Respiratory protection: Type approved RPE for dusts if required. Hand protection: rubber gloves. Eye protection: Goggles. Body protection: Overalls and boots		
9. Physical and chemical properties Appearance Form : Powder Colour : Pink Odour : Distinct Data relevant to safety Flash point : N/A Melting Point Range: 150 - 230 EC Ignition temperature N/A Density : 550 - 750 Kg / m3 Explosive Limits : N/A Solubility in water: Forms paste with water. Boiling Point / Range : N/A Vapour Pressure: N/A		

