



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 09/98/MW	Status	09/98/MW
Substance code:		Page :	1 (3)
1. Identification of the substance / preparation and company			
Product details			
Trade Name: DENLUBE (UNILUBE)			
Supplier details :			
PSP Dental Company Limited		Information provided by:	
3-5 Dylan Road		PSP Dental Company Limited	
Belvedere			
Kent DA17 5QS			
Emergency telephone number: 020-8311-7337			
Emergency fax number: 020 8310 0920			
2. Compostion/information on ingredients			
Chemical characterization			
Hazardous Components			
Ingredient	%	CAS No.	S Phrase R Phrase
Lubricant Oil	>50		
White Spirit	1.5		
Butane / Propane	20 - 50		
3. Hazards identification			
Unlikely to cause harmful effects under normal conditions of handling and use.			
Practice good hygiene. Wear gloves, overalls.			
4. First aid measures			
After contact with skin:			
Mild irritant, prolonged contact may cause defatting and dermatitis. Wash skin first with water and then, soap and water. If in doubt seek medical help.			
After contact with eyes:			
Mild irritant, wash open eye with of clean water for 10 minutes. Seek medical help.			
After ingestion:			
Seek medical advice.			
After inhalation:			
Remove to fresh air. If breathing is irregular or stopped, give artificial respiration. Give nothing by mouth. Seek medical help.			



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no.	Version 01/98/HLS	Status 01/98/HLS
Substance code: 10552		Page : 2 (3)
5. Fire-fighting measures: Fire exposed containers should be sprayed with water to lessen risk of explosion.		
6. Accidental release measures: Environmental precautions: Fire will cause dense black smoke, require self contained breathing apparatus. Do not run into drains. Methods for cleaning up / taking up: Absorb unto inert filler such as diatomaceous earth. Scoop into bag and seal.		
7. Handling and storage Handling: None required. Storage: Store in well ventilated cool area.		
8. Exposure controls / personal protection Hygiene measures Hand protection: rubber gloves. Eye protection: Goggles. Body protection: Air fed respiratory apparatus.		
9. Physical and chemical properties Appearance Form :Spray Colour :Clear Odour :Hydrocarbon Data relevant to safety Flash point : -40EC Ignition temperature: N/A Explosive Limits : >230EC auto flammable Boiling Point / Range : 154EC Melting Point Range: <0 Density : N/A Solubility in water: Insoluble in water. Vapour Pressure: N/A		

