



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no. Exactasil Putty	Version 09/98/MW	Status	09/98/MW	
Substance code:		Page :	1 (3)	
1. Identification of the substance / preparation and company				
Product details				
Trade Name: Exactasil Putty				
Supplier details :				
PSP Dental Company Limited		Information provided by:		
3-5 Dylan Road		PSP Dental Company Limited		
Belvedere				
Kent DA17 5QS				
		Emergency telephone number: 020-8311-7337		
		Emergency fax number: 020 8310 0920		
2. Composition/information on ingredients				
Chemical characterization				
Hazardous Components				
Ingredient	%	CAS No.	S Phrase	R Phrase
Talcum		14807-96-6		R96Xi
White Oil		8042-47-5		
3. Hazards identification				
Unlikely to cause harmful effects under normal conditions of handling and use.				
Practice good hygiene. Wear gloves, overalls.				
4. First aid measures				
After contact with skin:				
Mild irritant, prolonged contact may cause defatting and dermatitis. Wash skin first with water and then, soap and water. If in doubt seek medical help.				
After contact with eyes:				
Mild irritant, wash open eye with of clean water for 10 minutes. Seek medical help.				
After ingestion:				
Seek medical advice.				
After inhalation:				
Remove to fresh air. If breathing is irregular or stopped, give artificial respiration. Give nothing by mouth. Seek medical help.				



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no. Exactasil Putty Substance code: 10552	Version 01/98/HLS	Status 01/98/HLS Page : 2 (3)
5. Fire-fighting measures:		
6. Accidental release measures: Environmental precautions Prevent entry into drains and water courses. Avoid creating dusts. Methods for cleaning up / taking up: Wear approved respirator suitable for chemical dust concentrations encountered and impervious gloves (for bulk spillages). Sweep spill into containers for disposal. Wash contaminated area, wash affected skin areas with soap and water.		
7. Handling and storage Handling: Good General Ventilation Follow instructions enclosed in package and use only for purposes directed. Avoid contact with eyes. Avoid creating dusts. Storage: The container should be stored in normal dry place with lid firmly closed, since material is moisture sensitive.		
8. Exposure controls / personal protection Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work. Respiratory protection: Type approved RPE for dusts if required. Hand protection: rubber gloves. Eye protection: Goggles. Body protection: Overalls and boots		
9. Physical and chemical properties Appearance Form : Silicone Putty Colour : White Odour : slightly mint Data relevant to safety Flash point : Melting Point Range: N/A Ignition temperature : N/A Density : N/A Explosive Limits : N/A Solubility in water: Forms paste with water. Boiling Point / Range : N/A Vapour Pressure:		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no. Exactasil Putty	Version 01/98/HLS	Status 01/98/HLS
Substance code:		Page : 3 (3)
10. Stability and reactivity Thermal decomposition / Conditions to avoid Stable in normal conditions Hazardous reactions / Substances to avoid Unstable above 98EC gives Acrylate Vapour		
11. Toxicological information There is no Toxicological data available for this product. There is no known health risk if used properly. Prolonged skin and eye contact should be avoided. Remarks		
12. Ecological information Further ecological information Remarks		
13. Disposal considerations Product Substance via authorized waste disposal contractor to an approved waste disposal site, observing all local and national regulations.		
14. Transport information UN No: IMO: IATA: Correct Technical Name: ADR/RID: Not regulated for transport IMDG Class: Packaging Group: Packaging Group:		
15. Regulatory information: N/A Hazard symbols: R phrases: S phrases: S25: Avoid contact with eyes Other Regulations:		
16. Other information This Information is given in good faith, being based on the latest knowledge available to P.S.P.Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P.Dental Co. Ltd prior to use for further advice.		