



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no. Tray Tack Aerosol Substance code:	Version 02/98/TTA	Status 02/98/TTA Page : 1 (4)																																								
1. Identification of the substance / preparation and company Product details Trade Name: Tray Tack Aerosol Supplier details : PSP Dental Company Limited 3-5 Dylan Road Belvedere Kent DA17 5QS Information provided by: PSP Dental Company Limited Emergency telephone number: 0208- 311- 7337 Emergency fax number: 0208- 310- 0920																																										
2. Composition/information on ingredients Chemical characterization Hazardous Components <table border="1"><thead><tr><th>Ingredient</th><th>%</th><th>Cas No.</th><th>S Phrase</th><th>R Phrase</th></tr></thead><tbody><tr><td>Polyamide Resin</td><td>40</td><td></td><td>26 - 28</td><td>36/37</td></tr><tr><td>Xylene</td><td>25</td><td>1330-20-7</td><td>25</td><td>10-20/21-38</td></tr><tr><td>Isopropanol</td><td>35</td><td>67-63-0</td><td>7-16-24/25-26</td><td>11-36-67</td></tr></tbody></table> Non Hazardous Components <table border="1"><thead><tr><th>Ingredient</th><th>%</th><th>Cas No.</th><th>S Phrase</th><th>R Phrase</th></tr></thead><tbody><tr><td>Pigments</td><td><1</td><td></td><td></td><td></td></tr><tr><td>Butane</td><td></td><td>106-97-8</td><td>9-16</td><td>12</td></tr><tr><td>Propane</td><td></td><td>74-98-6</td><td>9-16</td><td>12</td></tr></tbody></table>			Ingredient	%	Cas No.	S Phrase	R Phrase	Polyamide Resin	40		26 - 28	36/37	Xylene	25	1330-20-7	25	10-20/21-38	Isopropanol	35	67-63-0	7-16-24/25-26	11-36-67	Ingredient	%	Cas No.	S Phrase	R Phrase	Pigments	<1				Butane		106-97-8	9-16	12	Propane		74-98-6	9-16	12
Ingredient	%	Cas No.	S Phrase	R Phrase																																						
Polyamide Resin	40		26 - 28	36/37																																						
Xylene	25	1330-20-7	25	10-20/21-38																																						
Isopropanol	35	67-63-0	7-16-24/25-26	11-36-67																																						
Ingredient	%	Cas No.	S Phrase	R Phrase																																						
Pigments	<1																																									
Butane		106-97-8	9-16	12																																						
Propane		74-98-6	9-16	12																																						
3. Hazards identification Irritating to eyes, slight skin irritation, Prolonged inhalation may cause dizziness.																																										
4. First aid measures after contact with skin Wash skin thoroughly with soap and water or proprietary cleaner. Do not use solvents or thinners. after contact with eyes Irrigate copiously for at least 10 minutes with clean fresh water. Seek medical attention. after ingestion Drink copious amounts of water. Seek immediate medical attention. Do not induce vomiting. after inhalation Remove to fresh air Keep patient warm and allow to rest.																																										



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no. Tray Tack Aerosol	Version 02/98/TTA	Status 02/98/TTA
Substance code:		Page : 2 (4)
5. Fire-fighting measures Suitable extinguishing media. Alcohol resistant foam, Carbon Dioxide, Powder, Water Spray/Mist. Do not use water jet.		
6. Accidental Release Measures. Environmental precautions Shut off sources of ignition. Ventilate area, Do not allow product to enter drains or water courses. Methods for cleaning up/taking up: Soak up with absorbent material (e.g. sand or earth) and place in a container for Disposal in accordance with local regulations		
7. Handling and storage Handling Avoid skin and eye contact. Avoid inhalation of vapour and spray mist. Do not use near sources of ignition. Storage Store at room temperature in a dry, well ventilated area away from sources of heat, Ignition and direct sunlight.		
8. Exposure controls / personal protection Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work. Respiratory protection :Unlikely to be necessary Hand protection : Use of protective gloves is recommended Eye protection :Wear goggles Body protection :Wear cotton or cotton/synthetic overalls		
9. Physical and chemical properties Appearance Form : Aerosol Colour : Red Odour : Aromatic Hydrocarbon Solvent Data relevant to safety Flash point 28EC Ignition temperature Explosive Limits 0.8-12.0 (% vol in air) Boiling Point / Range 82Ec Melting Point Range:<0Ec Density 0.86 Solubility in water Soluble Vapour Pressure		



Dental Company Limited

3/5 DYLAN ROAD
BELEVEDERE
KENT
DA17 5QS

HEALTH AND SAFETY DATA SHEET

Product no. Tray Tack Aerosol	Version 02/98/TTA	Status 02/98/TTA
Substance code:		Page : 3 (4)
10. Stability and reactivity Thermal decomposition/Conditions to avoid: Product is stable. Avoid naked flames, Red hot surfaces and other high temperature sources. Hazardous reactions/ Substances to avoid: Avoid oxidising agents, Strong Acids, Strong Alkalis.		
11. Toxicological information No Toxicological information available for this product. Remarks Avoid skin and eye contact.Avoid inhalation of vapours		
12. Ecological information No data available Further ecological information Remarks Do not allow product to enter drains and water courses		
13. Disposal considerations Product and empty containers should be disposed of in accordance with local and national regulations.		
14. Transport information Land transport / lorry / international / remarks ADR/RID 2 / 5F / 201 GGVS/GGVE 2 / 5F / 2201 aerosol HI-No UN-NO 1950 Inland waterways transport / remarks AND/ADNR: 2 5F aerosol AND/ADNR: 2 5F aerosol Marine transport / remarks IMDG/GGVSee-Code 2 EmS 2-13 MFAG 311 Marine pollutant - aerosols, limited quantities IMDG-PACK.GROUP: UN-Nr: 1950 Air transport / remarks ICAO/IATA-DGR: 2.1 UN no. : 1950 ICAO-PACK.GROUP: aerosols, flammable		

◆◆

Dental Company Limited

Product no. Tray Tack Aerosol Substance code:	Version 02/98/TTA	Status 02/98/TTA Page : 4 (4)
15. Regulatory information Hazard symbols R phrases : R 10-20/21-38 Highly flammable. Harmful by inhalation and in contact with skin. Irritating to skin. S phrases : S2-23-25-46 Keep out of reach of children. Do not breathe vapour/spray, Avoid contact with eyes. If Swallowed seek medical advice immediately and show container or label. Other regulations		
16. Other information This Information is given in good faith, being based on the latest knowledge available to P.S.P .Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P .Dental Co. Ltd prior to use for further advice.		