



Dental Company Limited

3/5 DYLAN ROAD  
BELEVEDERE  
KENT  
DA17 5QS

### HEALTH AND SAFETY DATA SHEET

Product no. Tray Tack Liquid	Version 09/98/MW	Status 02/98/DRG
Substance code:	Page :	1 (3)
<b>1. Identification of the substance / preparation and company</b>		
Product details		
Trade Name: Tray Tack		
Supplier details :		
PSP Dental Company Limited 3-5 Dylan Road Belvedere Kent DA17 5QS		Information provided by: PSP Dental Company Limited
<b>Emergency telephone number: 020-8311-7337</b> <b>Emergency fax number: 020 8310 0920</b>		
<b>2. Compostion/information on ingredients</b>		
Chemical characterization		
Hazardous Components		
Ingredient	%	CAS No. S Phrase R Phrase
Polyamide Resin	40	
Xylene	25	1330-20-7
Isopropanol	35	67-63-0
Non Hazardous Components		
Ingredient	%	CAS No
Pigments	<1	
<b>3. Hazards identification</b>		
<b>Irritating to eyes, slight skin irritant, prolonged inhalation may cause dizziness</b>		
<b>4. First aid measures</b>		
After contact with skin: Wash skin with water and soap or proprietary cleaner . Do not use solvents or thinners.		
After contact with eyes: Irrigate copiously for at least 10 minutes with clean fresh water. Seek medical help.		
After ingestion: Drink copious amounts of water. Seek immediate medical attention.		
After inhalation: Remove to fresh air. Keep patient warm and allow to rest.		



Dental Company Limited

3/5 DYLAN ROAD  
BELEVEDERE  
KENT  
DA17 5QS

### HEALTH AND SAFETY DATA SHEET

Product no. Tray Tack Liquid	Version 01/98/HLS	Status 01/98/HLS
Substance code:	Page : 2 (3)	
<b>5. Fire-fighting measures:</b> Suitable extinguishing media. Alcohol resistant foam, Carbon Dioxide Powder, Water Spray / mist Do not use water jet.		
<b>6. Accidental release measures:</b> Environmental precautions Shut off sources of ignition. Ventilate area. Do not allow product to enter drains or water courses.  Methods for cleaning up / taking up: Soak up with absorbent material (e.g sand or earth) and place in a container for disposal in accordance with local Regulations.		
<b>7. Handling and storage</b> Handling: Avoid skin and eye contact. Avoid inhalation of vapour and spray mist. Do not use near sources of ignition.  Storage: Store at room temperature in a dry, well ventilated area away from sources of heat, ignition and direct sunlight.		
<b>8. Exposure controls / personal protection</b>  Hygiene measures Remove soiled or soaked clothing immediately. Wash hands before breaks and after work.  Respiratory protection: Unlikely to be necessary  Hand protection: Use of protective gloves is recommended  Eye protection: Wear Goggles.  Body protection: Wear cotton or cotton / synthetic overalls		
<b>9. Physical and chemical properties</b> <b>Appearance</b> Form : Liquid Colour : Red Odour : Aromatic Hydrocarbon Solvent <b>Data relevant to safety</b> Flash point : 28°C Melting Point Range: <0 °C  Ignition temperature N/A Density : 0.86  Explosive Limits : 1.1 – 12.0 (% vol in air) Solubility in water: Insoluble  Boiling Point / Range : 82°C Vapour Pressure:		



Dental Company Limited

3/5 DYLAN ROAD  
BELEVEDERE  
KENT  
DA17 5QS

### HEALTH AND SAFETY DATA SHEET

Product no. Tray Tack Liquid	Version 01/98/HLS	Status 01/98/HLS
Substance code:	Page :	3 (3)
<b>10. Stability and reactivity</b> Thermal decomposition / Conditions to avoid Avoid naked flames, red hot surfaces and other high temperature sources  Hazardous reactions / Substances to avoid Avoid oxidising agents, strong acids, strong alkalis.		
<b>11. Toxicological information</b> There is no Toxicological data available for this product. There is no known health risk if used properly. Prolonged skin and eye contact should be avoided.  Remarks Avoid skin and eye contact. Avoid inhalation of vapours		
<b>12. Ecological information</b> Further ecological information Remarks Do not allow product to enter drains and water courses		
<b>13. Disposal considerations</b> Product and empty containers should be disposed of in accordance with local and national regulations.		
<b>14. Transport information</b>  UN No: 1993 IMO: IATA:  IMDG Class: 3230 Packaging Group: Packaging Group:  Correct Technical Name: ADR/RID:		
<b>15. Regulatory information</b> Hazard symbols: R phrases: R11-20/21-38. Highly flammable. Harmful by inhalation and in contact with skin. Irritating to skin  S phrases: S2-23-25-46. Keep out of reach of children. Do not breathe vapour / spray. Avoid contact with eyes. If swallowed seek medical advice immediately and show this container or label.  Other Regulations:		
<b>16. Other information</b> This Information is given in good faith, being based on the latest knowledge available to P.S.P.Dental Co. Ltd. No known relevant information has been omitted from the product information, risk information and safety precaution sheets. The user should not assume therefore, the product is suitable for any abnormal use. If the information provided is insufficient to ensure safety in any particular application, please contact P.S.P.Dental Co. Ltd prior to use for further advice.		